


Dana 2.2 MG G

Automatic Biopsy gun for disposable Tru-Cut type biopsy needles

* Real size



- Automatic Biopsy gun for needles Tru-Cut type with 0mm or 22mm penetration.
- Double trigger, frontal and back.
- Compatible with Histo Biocore II MG biopsy needles. 220 gr. Weight.
- Visual and tactile safety.
- Reusable system for harvesting soft tissue specimens.
- Aluminium structure.
- Easy sample retrieval without removing needle from the gun.
- Easily opened.
- It allows sampling in axillary lymph nodes; Breasts with prosthesis; Zones of risk with no place for shots with advance.

| CODE | COLOR | DESCRIPTION |
|--------|---|--|
| 360015 |  | DANA 2.2 MGG, penetration depth 0 mm 22 mm |

* This system can use needles from, 12, 14, 16 y 18 Ga.

"Certified company ISO 13485:2016 / ISO 9001:2015"



histo

Av. General Paz 315 -Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



GESTION
DE LA CALIDAD
ISO 13485:2016
ISO 9001:2015



anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica

Biocore II MG

Core biopsy needle for soft tissue specimen retrieval. TRU-CUT type compatible with various automatic biopsy guns commercially available

- Core biopsy needle for soft tissue specimen retrieval. TRU - CUT type compatible with Dana 2.2 MG and various instruments commercially available.
- Colored coded hubs to identify needle gauge.
- Centimeter depth markings for greater accuracy.
- Exclusive "spacer" system to introduce needle in the lesion before insertion into the biopsy instrument (gun).
- "Spacer" will detach after placing the needle into the instrument (gun).
- Hubs adapted to allow syringe use.



| CODE | COLOR | GAUGE | LENGTH |
|--|-------|----------------------------|---|
| 551210 551213 | ● | 12 12 | 10 cm 13 cm |
| 551410 551413 551415 551420 | ● | 14 14 14 14 | 10 cm 13 cm 15 cm 20 cm |
| 551610 551613 551615 551620 551625 | ● | 16 16 16 16 16 | 10 cm 13 cm 15 cm 20 cm 25 cm |
| 551810 551815 551820 551825 551830 | ● | 18 18 18 18 18 | 10 cm 15 cm 20 cm 25 cm 30 cm |

* Disposable product, for single patient use only.

* See instructions contained in package needle for correct use.

"Certified company ISO 13485:2016 / ISO9001:2015 "

HISTO

Av. General Paz 315 -Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



GESTION
DE LA CALIDAD
ISO 13485:2016
ISO 9001:2015



anmat
Administración Nacional de Medicamentos,
Biotecnología y Tecnología Médica



Biocore II BR

Core biopsy needle for soft tissue specimen retrieval. TRU-CUT type compatible with various automatic biopsy guns commercially available

- Core biopsy needle for soft tissue specimen retrieval. TRU - CUT type compatible with various instruments commercially available.
- Colored coded hubs to identify needle gauge.
- Centimeter depth markings for greater accuracy.
- Exclusive “spacer” system to introduce needle in the lesion before insertion into the biopsy instrument (gun).
- “Spacer” will detach after placing the needle into the instrument (gun).
- Hubs adapted to allow syringe use.

| CODE | COLOR | GAUGE | LENGTH |
|--|-------|----------------------------|---|
| 661410 661415 661420 | ● | 14 14 14 | 10 cm 15 cm 20 cm |
| 661610 661615 661620 661625 | ● | 16 16 16 16 | 10 cm 15 cm 20 cm 25 cm |
| 661810 661815 661820 661825 661830 | ● | 18 18 18 18 18 | 10 cm 15 cm 20 cm 25 cm 30 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.

“Certified company ISO 13485:2003 / ISO 9001:2008”

HISTO

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar

anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica



Biocore II

Core biopsy needle for soft tissue specimen retrieval. TRU-CUT type compatible with various automatic biopsy guns commercially available

- Core biopsy needle for soft tissue specimen retrieval. TRU - CUT type compatible with various instruments commercially available.
- Colored coded hubs to identify needle gauge.
- Centimeter depth markings for greater accuracy.
- Hubs adapted to allow syringe use.
- 19 mm Notch.

| CODE | COLOR | GAUGE | LENGTH |
|--------|-------|-------|--------|
| 441410 | ● | 14 | 10 cm |
| 441415 | | 14 | 15 cm |
| 441420 | | 14 | 20 cm |
| 441610 | ● | 16 | 10 cm |
| 441615 | | 16 | 15 cm |
| 441620 | | 16 | 20 cm |
| 441625 | | 16 | 25 cm |
| 441810 | ● | 18 | 10 cm |
| 441815 | | 18 | 15 cm |
| 441820 | | 18 | 20 cm |
| 441825 | | 18 | 25 cm |
| 441830 | | 18 | 30 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.

"Certified company ISO 13485:2003 / ISO 9001:2008"

HISTO

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
 Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar

anmat
 Asociación Nacional de Fabricantes de Instrumentos Médicos



Biocore II TC

Aguja de biopsia semiautomática
tipo Trucut

- Operable con una mano.
- Dos opciones de carga y disparo.
- Desarmable modo coaxial.
- 10 o 20 mm de penetración.
- Calibres disponibles: 12GA; 14GA; 16GA; 18GA.
- Código de producto " 77 XX XX "



Escaneá el código
QR para obtener
más información
del producto.



* Producto descartable, para utilizar una única vez.

* Para su uso correcto, ver información en el instructivo de uso junto al embalaje del producto.

HISTO

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



GESTIÓN
DE LA CALIDAD
ISO 13485:2016
ISO 9001:2015



anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica



Blo-K

Prostate pre biopsy blocking needle

- Used through the ultrasound transducer prior to inserting the needle biopsy.
- Allows you to use all kinds of syringes.
- The bezel has a special design for chuck not need to be introduced.
- It has depth markings every inch.
- Terminal has a translucent plastic, for better viewing.

| CODE | COLOR | GAUGE | LENGTH |
|--------|-------|-------|--------|
| 812020 | ○ | 20 | 20 cm |
| 812025 | | 20 | 25 cm |

* Disposable product, for single patient use only.

* See instructions contained in package needle for correct use.

"Certified company ISO 13485:2003 / ISO 9001:2008"

HISTO

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar

anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica



Mamawire Plus

Breast localization needles

- Breast localization needle “double fixing effect”.
- Provides maximum retention avoiding wire migration.
- Specially recommended for large breast or “loose” breast tissue.
- Ultrasound and x-ray detectable.
- Introducer needle depth marks allows to calculate lesion distance penetration.
- Localization wire has a mark to signify the point at which the barb will be deployed.
- End hole needle design allows dye injection or fluid aspiration.

| CODE | INTRODUCER NEEDLE | LOCALIZATION WIRE |
|--------|-------------------|-------------------|
| 312005 | 20 Ga x 5 cm | 15 cm |
| 312007 | 20 Ga x 7 cm | 17 cm |
| 312010 | 20 Ga x 10 cm | 20 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.

“Certified company ISO 13485:2003 / ISO 9001:2008”

histo

Av. General Paz 315 -Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



Mamawire

Breast localization needles



- Breast localization needle.
- Ultrasound and x-ray detectable.
- Introducer needle depth marks allows to calculate lesion distance penetration.
- Localization wire has a mark to signify the point at which the barb will be deployed.
- End hole needle design allows dye injection or fluid aspiration.

| CODE | INTRODUCER NEEDLE | LOCALIZATION WIRE |
|--------|-------------------|-------------------|
| 212005 | 20 Ga x 5 cm | 15 cm |
| 212007 | 20 Ga x 7 cm | 17 cm |
| 212010 | 20 Ga x 10 cm | 20 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.

"Certified company ISO 13485:2003 / ISO 9001:2008"

histo

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar

anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica





Mamawire Fix D

Landmark for intervened lesions

- After biopsy with vacuum method or Core, may, through the left channel made by the needle biopsy, enter the MAMAWIRE FIX and reach the desired depth pressing the actuator and leaving the clip attached.
- Without performing a biopsy allows the MAMAWIRE FIX to be introduced, after reaching the desired depth, push the actuator depositing the clip.
- Ultrasound, X-ray, and nuclear magnetic resonance detectable.
- Depth markings every inch, and has a sure shot for the safe handling of ejecting the clip.
- Clip reduced very easy and quick application size.
- Shows the area where the histological biopsy was performed for later tissue controls.
- Tip with patient access capability.
- Indicates breast lesions without biopsy for future controls.



| CODE | COLOR | GUAGE | LENGTH |
|--------|-------|-------|--------|
| 911810 | ○ | 18 | 10 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.

"Certified company ISO 13485:2016 / ISO 9001:2015"

HISTO

Av. General Paz 315 - Villa Raffo, Bs. As. - AR - 1674
 Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



GESTION
DE LA CALIDAD
ISO 13485:2016
ISO 9001:2015



anmat
Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica





Mamawire Fix

Landmark for intervened lesions

- After biopsy with vacuum method or Core, may, through the left channel made by the needle biopsy, enter the MAMAWIRE FIX and reach the desired depth pressing the actuator and leaving the clip attached.
- Ultrasound, X-ray, and nuclear magnetic resonance detectable.
- Depth markings every inch, and has a sure shot for the safe handling of ejecting the clip.
- Clip reduced very easy and quick application size.
- Shows the area where the histological biopsy was performed for later tissue controls.



| CODE | COLOR | GUAGE | LENGTH |
|--------|-------|-------|--------|
| 911810 | ○ | 18 | 10 cm |

- * Disposable product, for single patient use only.
- * See instructions contained in package needle for correct use.t

"Certified company ISO 13485:2016 / ISO9001:2015 "

HISTO

Av. General Paz 315 -Villa Raffo, Bs. As. - AR - 1674
 Tel./Fax: (005411) 4712 5941 // 4757 1437
www.histo.com.ar // info@histo.com.ar



GESTION DE LA CALIDAD
 ISO 13485:2016
 ISO 9001:2015



anmat
 Administración Nacional de Medicamentos,
 Alimentos y Tecnología Médica

